



# Oxygen Sensor


## Product description

Zimmer and Peacock's Oxygen sensors can be used in continuous or discrete mode. The ZP Oxygen sensors are intended for anyone looking for a robust low-cost platform for oxygen sensing.



Scan the QR-code for videos on how to use ZP sensors!

## Performance



ZP recommends the Anapot potentiostat when using our sensors.

The software "Lisa Analytics" for the Anapot can be found at:

<https://www.zimmerpeacocktech.com/knowledge-base/software/>

The link also contains videos on how to properly configure your device with the software.

## Sensor materials


Substrate	Reference electrode	Working electrode	Counter electrode
Alumina	Silver/Silver Chloride (60/40)	Platinum	Silver/Silver Chloride (60/40)

[www.zimmerpeacock.com](http://www.zimmerpeacock.com)



# Oxygen Sensor


## Mechanical drawings



Thickness = 0.625

All units in mm

Sensor product code: A-AC-PAC-303-O



## Storage

Storage: Room temperature 20 - 50% RH

Protect from exposure to UV-light

Recommended usage within 6 months of fabrication

## Safety

Take caution when handling the sensors as the edges and corners are sharp

ZP recommends using gloves

This product is not tested for biocompatibility and ZP takes no responsibility for in-vivo usage

## Customizability

Zimmer and Peacock can also make customized sensors with the option to target other analytes than those listed in this datasheet. We can offer different electrode configurations, geometry, and materials. Please contact us through the contact form on [www.zimmerpeacock.com](http://www.zimmerpeacock.com) or by e-mail on [sales@zimmerpeacock.com](mailto:sales@zimmerpeacock.com) for questions regarding customized sensors

[www.zimmerpeacock.com](http://www.zimmerpeacock.com)